

Our Sights Towards a Greener Tomorrow...

...Taking pride while being responsible for producing environmentally safe products, that will enable us to achieve our goals towards a *Healthier and Greener future...*

ECO-BASE™ SURGICAL INSTRUMENT DETERGENT

Eco-Base Surgical Instrument Detergent is a non-enzymatic detergent designed for automated processes, such as tunnel washers, ultrasonic washers and cart washers. It is an effective detergent used for the cleaning of metallic and plastic medical, surgical and dental instruments, including flexible and rigid endoscopes.

What are the Benefits?

Effective & Convenient: Eco-Base Surgical Instrument Detergent is able to effectively clean surgical instrumentation in all automatic equipment.

Performs Well in Hard Water: Formulated to perform well in hard water.

PT.# 11050-04 : 4/GALLONS/CASE

PT.# 11050-2.5 : 2/2.5 GALLONS/CASE

PT.# 11050-05 : 5/GALLON DRUM

Low Foaming: Specifically designed for automatic machines, such as ultrasonic cleaners and endoscope reprocessors, as per manufacturer's recommendations.

Neutral pH: pH is neutral once diluted.

PT.# 11050-15 : 15 GALLON DRUM

PT.# 11050-30 : 30 GALLON DRUM



ECOZYME™ ULTRA CONCENTRATE SPRAY

Ecozyme Spray is an effective multi-enzymatic spray used immediately after surgery to begin the breakdown of organic and inorganic biological residue on surgical instrumentation. Ecozyme Spray avoids the use of water in the O.R. It is used for the cleaning of metallic and plastic medical, surgical and dental instruments, including flexible and rigid endoscopes.

What are the Benefits?

Fast Acting: Prevents organic biological residue from drying onto surgical instruments

Convenient: Ecozyme Spray is compatible with medical, surgical, and dental instruments. Rigid and Flexible ENDOSCOPES.

Effective: Works on contact

PT.# 27545-16 : 12/16OZ CONTAINERS/CASE

PT.# 27545-32 : 6/32OZ CONTAINERS/CASE

PRY® SUPER STRENGTH CLEANING CREAM

Pry Instrument Cleaning Cream is an all purpose cleaner for the removal of stains, discoloration, and rust from metal instrumentation, sterilizer doors and internal washer chambers.

What are the Benefits?

Effective & Convenient: Pry Instrument Cleaning Cream is able to effectively remove stains, discoloration, and rust from metal instrumentation, sterilizer doors and internal washer chambers.

PT.# 43394-91



* All Detergents Are Classified As "Ultimate Biodegradable" OCED Testing Standards

Visit us at: www.prolinesolutions.com

Email us at: prolinesolutions@aol.com



**PRACTICE
Greenhealth**

Meets Practice Greenhealth Specifications

